

MegaMover® Portable Transport Chair



NEW
Yellow handles
easily identify
Transport Chair



The MegaMover® Transport Chair is a strong, disposable alternative for transporting patients down stairs and through narrow hallways or extracting them from confined spaces like bathrooms.

- **Sanitary and Disposable.** If it gets soiled just throw it away.
- Made of a fluid-resistant material strengthened with reinforced nylon straps.
- **Strong.** Carries up to 350 pounds with the weight distributed across eight nylon reinforced handles.
- Carrying handles give a firm grip for up to four people to share the load.
- **Compact.** Weighs only 10 ounces and rolls to the size of a water bottle.



From the makers of
MEGA MOVER®
portable transport unit

Product Order Information

Code	Product Name	Material	Color	Size	Packing
56307	MegaMover® Transport Chair	Nonwoven/Poly	White/Yellow	35" x 35"	10/Case
53376	MegaMover® Plus	Nonwoven/Poly	Gray/Black	40" x 80"	10/Case
51926	MegaMover®	Nonwoven/Poly	White/Black	40" x 80"	10/Case