

# MEGAMOVER® TACTICAL

PORTABLE TRANSPORT UNIT



## MEGAMOVER® TACTICAL

- Designed to be dragged or carried--up to 500 lbs.
- Very lightweight, ultra compact, strong and rugged.
- Preferred by small Rescue Task Force Teams.
- Requires only one to two responders per victim.
- Can be loaded from either direction.
- Center strap with buckle to secure victim; drag straps on both ends.
- Patent pending.
- Visit [grahammedical.com](http://grahammedical.com) for other MegaMover® products.



Item #	Material	Color	Size	Pack
81644	Nonwoven/Poly	Gray	24" x 48"	10/CS






# MEGAMOVER® DISASTER RESPONSE

PORTABLE TRANSPORT UNIT



## MEGAMOVER® DISASTER RESPONSE

- Designed to be carried or dragged--up to 500 lbs.
- Lightweight, compact, strong and rugged.
- Used by Rescue Task Force Teams with two to eight responders per victim.
- Screen printed with "Head" and "Foot" to improve efficiency under pressure.
- Under arm strap with buckle offset to secure victims.
- Patent pending.
- Visit [grahammedical.com](http://grahammedical.com) for other MegaMover® products.



Item #	Material	Color	Size	Pack
81707	Nonwoven/Poly	Gray	32" x 72"	8/CS

