

# Snap-Closure Gown

## Disposable Snap-Closure Gown Ideal For Accessibility

- Easy access snap-closure sleeves and back.
- Made of soft, linen-like material that is fluid-resistant, static-free, lightweight and breathable.
- Greatly reduce the risk of cross-contamination with a new, clean gown every time.
- Safe to use in X-Ray, MRI and Autoclave machines.



| Code  | Item | Description                | Material | Color  | Trim  | Size  | Size Info | Pack  |
|-------|------|----------------------------|----------|--------|-------|-------|-----------|-------|
| 84792 | Gown | Snap Sleeves, Back-Opening | Nonwoven | Maroon | White | S     | 20" x 34" | 30/CS |
| 84793 | Gown | Snap Sleeves, Back-Opening | Nonwoven | Maroon | White | L     | 30" x 46" | 30/CS |
| 84794 | Gown | Snap Sleeves, Back-Opening | Nonwoven | Maroon | White | 4/5XL | 42" x 50" | 30/CS |