



# ***EMS DISPOSABLES***

*DURABLE*

*RESPONSIVE*

*VERSATILE*




graham  
medical<sup>®</sup>  
The Promise of Protection<sup>®</sup>

# FITTED *STRETCHER SHEETS & EQUIPMENT COVERS*

Multiple sizes and fits to accommodate a variety of applications.



SnugFit® and PowerFit® sheets have elastic corners to prevent slippage when patient moves.



The SnugFit® Equipment Cover #72930 is custom designed to fit the Stryker® XPS Mattress.



**Graham Medical's fitted sheets and equipment covers are quick and easy to use, and fit almost all equipment.**

- Reduce the risk of infection with fluid-resistant and fluid-barrier designs.
- Covers durable enough to aid patient transfer up to 400 lbs.
- Keep equipment contamination free and minimize cleaning time.
- Latex-free.

### **PocketFit®**

- No elastic allows sheet to slip on and off mattresses quickly and easily.
- Fluid-resistant.

### **SnugFit®**

- The elastic, fitted design keeps sheet tight and flat on stretcher to prevent slippage during patient movement.
- Easy, secure fit for standard stretchers and gurneys.
- Available in fluid-resistant and barrier.

### **PowerFit®**

- Same great features as SnugFit®.
- Fits oversized stretchers and beds.
- Available in fluid-resistant and barrier.



Code	Item	Material	Color	Size	Case	Transfer	Quality
65233	PocketFit® Equipment Cover	Nonwoven	Blue	30" x 72"	50/CS	200 lbs.	Standard
65245	PocketFit® Equipment Cover	Nonwoven	White	30" x 72"	50/CS	250 lbs.	Premium
65247	SnugFit® Equipment Cover	Nonwoven	Blue	30" x 84"	50/CS	NA	Standard
44547	SnugFit® Fitted Sheet	Nonwoven	Blue	30" x 84"	5/BG; 10BG/CS	250 lbs.	Standard
53378	SnugFit® Fitted Sheet - Barrier	Nonwoven/Poly	Blue	30" x 84"	5/BG; 10BG/CS	250 lbs.	Premium
49896	SnugFit® Fitted Sheet	Nonwoven	Blue	40" x 89"	5/BG; 5BG/CS	250 lbs.	Standard
54311	SnugFit® Fitted Sheet - Barrier	Nonwoven/Poly	Blue	40" x 89"	5/BG; 5BG/CS	250 lbs.	Premium
61234	SnugFit® Fitted Sheet - ED/ER Mattress	Nonwoven	Blue	40" x 89"	25/CS	250 lbs.	Premium
72930	SnugFit® Equipment Cover - Fits Stryker® XPS	Nonwoven	Blue	36" x 84"	5/BG; 6BG/CS	300 lbs.	Premium
54419	PowerFit® Fitted Sheet	Nonwoven	Blue	33" x 89"	5/BG; 10BG/CS	250 lbs.	Standard
54420	PowerFit® Fitted Sheet - Barrier	Nonwoven/Poly	Blue	33" x 89"	5/BG; 10BG/CS	250 lbs.	Premium
85910	PowerFit® Equipment Cover - Barrier Heavyweight	Nonwoven/Poly	White	33" x 89"	5/BG; 10BG/CS	400 lbs.	Premium

# FLAT SHEETS

Single-use flat sheets offer convenience for providers and protection for patients.

- Multi-layered tissue sheets provide maximum absorbency and comfort.
- Tissue/Poly sheets act as a powerful protective barrier.
- Nonwoven sheets bring the strength and comfort of linens to a convenient, disposable sheet.

	Code	Item	Material	Color	Size	Case	Quality
USA	34682	Flat Sheet	2-Ply Tissue	White	40" x 72"	50/CS	Economy
USA	304	Flat Sheet	3-Ply Tissue	White	40" x 72"	50/CS	Standard
USA	305	Flat Sheet	3-Ply Tissue	White	40" x 90"	50/CS	Standard
USA	52166	Flat Sheet	Tissue/Poly	White/Blue	40" x 90"	50/CS	Economy
USA	322	Flat Sheet	Tissue/Poly-Hvywt.	White/Blue	40" x 72"	50/CS	Standard
USA	323	Flat Sheet	Tissue/Poly-Hvywt.	White/Blue	40" x 90"	50/CS	Standard
USA	324	Flat Sheet - Fanfold	Tissue/Poly/Tissue	Blue	40" x 72"	48/CS	Standard
USA	329	Flat Sheet - Fanfold	Tissue/Poly/Tissue	Blue	40" x 84"	48/CS	Standard
USA	332	Flat Sheet - Fanfold	Tissue/Poly/Tissue	Blue	40" x 96"	48/CS	Standard
USA	334	Flat Sheet - Fanfold	Tissue/Poly/Tissue-Hvywt.	Blue	40" x 84"	48/CS	Standard
	65232	FlexDrape®	Nonwoven	Navy	40" x 84"	50/CS	Premium
USA	50588	FlexDrape® - Fanfold	Nonwoven	Blue	50" x 84"	50/CS	Premium



# PILLOWS & PILLOWCASES

Disposable pillows and pillowcases offer comfort for patients and versatility for providers.

## FlexAir® Pillows

- Can be used as a headrest and for limb or posture support.
- Space-saving design stores flat.
- Adjustable valve allows easy inflation or deflation in seconds.

## Pillowcases

- Tissue/Poly provides high absorbency and acts as a true fluid barrier.
- Nonwoven provides the same fluid barrier with optimal linen-like comfort.



	Code	Item	Material	Color	Size	Case	Quality
	50349	FlexAir® Inflatable Pillow	Nonwoven/Poly	White	14.5" x 10.5"	50/CS	Standard
	51038	FlexAir® Plus Inflatable Pillow	Nonwoven/Poly	White	19" x 12.5"	50/CS	Standard

	Code	Item	Material	Color	Size	Case	Quality
USA	48766	Pillowcase	Tissue/Poly	Blue	21" x 30"	100/CS	Economy
USA	47256	Pillowcase	Tissue/Poly	White	21" x 30"	100/CS	Economy
	53157	Solacel™ Pillowcase	Tissue/Extruded Poly	White	22" x 30"	100/CS	Standard
USA	354	Pillowcase	Nonwoven	White	21" x 27.75"	100/CS	Premium

USA = Made In The U.S.A.

# BLANKETS *SINGLE-USE*

Graham Medical's versatile blanket line provides a solution for every scenario.

- Reduce cross-contamination with single-use products.
- Sized to provide total coverage and protection.
- Comfort1® Blankets can be used in warmers up to 150 degrees.
- VisiBlankets® draw immediate attention to priority patients.
- Hypothermia blankets are a lightweight material with superior heat retention.

	Code	Item	Material	Color	Size	Case	Quality
USA	54849	VisiBlanket®	2-Ply Tissue/Poly	Yellow/White	54" x 84"	5/BG; 5BG/CS	Economy
USA	77700	VisiBlanket®	2-Ply Tissue/Poly	Yellow/White	56" x 90"	1/BG; 24BG/CS	Standard
	53382	VisiBlanket®	Poly/Airlaid	Yellow	50" x 84"	1/BG; 25BG/CS	Premium
	53377	Hypothermia Blanket	Mylar	Silver	52" x 84"	1/BG; 50BG/CS	Standard
	78313	Comfort1® Blanket	Polyester	Green	36" x 60"	1/BG; 10BG/CS	Premium
USA	84869	Comfort1® Blanket	Polyester	Gray	40" x 80"	1/BG; 25BG/CS	Economy
USA	62292	Comfort1® Blanket	Polyester	Gray	50" x 84"	1/BG; 10BG/CS	Standard
USA	52038	Comfort1® Blanket	Polyester	White	50" x 84"	1/BG; 10BG/CS	Premium
	77723	Comfort1® Polar Fleece	Fleece	Navy	58" x 80"	1/BG; 10BG/CS	Premium



# INSTAKIT® *BEDDING & TRANSPORT KITS*

Prepackaged bedding and transport kits save time when every second matters.

- Everything you need in one ready-to-go package.
- Eliminate the hassle of replacing ripped and stained linens.
- Save on storage space.

### 50983 InstaKit® Bedding Kit 25/CS

Components:	44547	SnugFit® Fitted Sheet 30"x84" Nonwoven Blue
	329	Flat Sheet 40"x84" Tissue/Poly/Tissue Blue
	360	Pillowcase 21"x30" Tissue/Poly White
	50984	Absorbent Underpad 17"x24" Blue

### USA 65381 InstaKit® Bedding Kit 25/CS

Components:	323	Flat Sheet 40"x90" Tissue/Poly-Heavyweight White/Blue
	324	Flat Sheet 40"x90" Tissue/Poly/Tissue-Heavyweight Blue
	47256	Pillowcase 21"x30" Tissue/Poly White

### 84600 InstaKit® Transport Kit 10/CS

Components:	77723	Comfort1® Polar Fleece Blanket 50"x84" Navy
	52037	MegaMover® 350i Heavy-Duty Barrier Sheet 20"x84" Gray



**CUSTOM KITS  
AVAILABLE**

Dependent On  
Volume & Components

USA = Made In The U.S.A.

# MEGAMOVER® PORTABLE TRANSPORT UNITS

**MegaMovers® are compact, cost-effective portable transport units used to transport, transfer or rescue patients from areas inaccessible to stretchers. The nonwoven/poly, latex-free construction is fluid resistant and folds or rolls up for easy storage. Some models are equipped with PowerGrip® handles for added comfort and durability.**

## MegaMover® – The Original

- Transport up to 1,000 lbs.
- 14 handles.
- Available with standard or PowerGrips®.

## MegaMover® Plus

- Backboard pockets.
- Transport up to 1,000 lbs.
- 14 handles.
- Available with standard or PowerGrips®.

## MegaMover® Select

- Transport up to 1000 lbs.
- 12 handles with multi-grip options.
- Two drag/pull straps on each end.

## MegaMover® 350i Heavy-Duty Barrier Sheet

- Fluid barrier sheet.
- Transport up to 350 lbs.
- Use under patient to protect equipment.
- Use over patient for cover.

## MegaMover® Transport Chair

- Ideal for use in confined space and when patients need to be in a seated position.
- Transport up to 350 lbs.
- 8 handles.
- Weighs only 10 oz.
- Compact, rolls to the size of a water bottle.

## SPECIALIZED MEGAMOVERS®

### MegaMover® Tactical

- Ideal for Tactical Active Shooter events.
- Weighs 1.3 lbs., easy to carry.
- Transport up to 500 lbs.
- 1 or 2 responders.
- Load patient from either direction.
- Secure patient with buckle strap.
- Drag from either end.

### MegaMover® Disaster Response

- Ideal for Disaster Evacuation events.
- Weighs 2 lbs., easy to carry.
- Transport up to 500 lbs.
- 10 handles with multi-grip options.
- 2 to 8 responders.
- Printed with “Head” and “Foot” guides to improve efficiency under pressure.
- Secure patient with underarm strap.



The MegaMover® Transport Chair is ideal for carrying patients in areas of confined space.




Also available with PowerGrips® for added comfort and durability.



MegaMover® Plus accommodates standard backboards for immobilization. *Backboard not included.*

## MEGAMOVER® portable transport unit

Code	Item	Material	Color	Webbing	Size	Case	Max Wt.	Quality
51926	MegaMover®	Nonwoven/Poly	White	Black	40" x 80"	10/CS	1000 lbs.	Standard
66563	MegaMover® With PowerGrips®	Nonwoven/Poly	White	Black	40" x 80"	8/CS	1000 lbs.	Premium
53376	MegaMover® Plus (Backboard Pockets)	Nonwoven/Poly	Gray	Black	40" x 80"	10/CS	1000 lbs.	Premium
65459	MegaMover® Plus With PowerGrips® (Backboard Pockets)	Nonwoven/Poly	Gray	Black	40" x 80"	8/CS	1000 lbs.	Premium
72657	MegaMover® Select	Nonwoven/Poly	White	Black/Green	40" x 80"	8/CS	1000 lbs.	Premium
52037	MegaMover® 350i Heavy-Duty Barrier Sheet	Nonwoven/Poly	Gray	NA	50" x 84"	25/CS	350 lbs.	Economy
56307	MegaMover® Transport Chair	Nonwoven/Poly	White	Yellow	36" x 36"	10/CS	350 lbs.	Standard
81644	MegaMover® Tactical	Nonwoven/Poly	Gray	Black	24" x 48"	10/CS	500 lbs.	Premium
81707	MegaMover® Disaster Response	Nonwoven/Poly	Gray	Black	32" x 72"	8/CS	500 lbs.	Premium





Graham Medical is the leader in the manufacturing and marketing of disposable EMS products.

For more information on Graham Medical and our products  
call Customer Service at **800-558-6765**,  
visit [grahammedical.com](http://grahammedical.com) or email [info@grahammedical.com](mailto:info@grahammedical.com)



2273 Larsen Road | P.O. Box 19031 | Green Bay, WI 54307-9031 | 800.558.6765 | [info@grahammedical.com](mailto:info@grahammedical.com) | [www.grahammedical.com](http://www.grahammedical.com)

©2020 Little Rapids Corporation



PM293E0220