

Procedure Masks

Our USA made masks are comfortable to wear with a **skin friendly contact layer**, ample coverage, and **soft, flexible ear loops**. Crafted of quality USA made materials, our single-use masks are **comfortable to wear** for long periods of time.

Metal Free Procedure Masks Series 1



- ✓ Contains no metal - safe for use in MRI Machines and Behavioral Health settings
- ✓ Soft, comfortable ear loops
- ✓ 3-Ply Construction
- ✓ Latex-free
- ✓ Skin-friendly contact layer that meets ISO Standard 10993 for skin irritation and sensitization



Procedure Masks Series 1 & 2



- ✓ 3-Ply Construction
- ✓ Soft, comfortable ear loops
- ✓ Adjustable nose piece for better fit
- ✓ Latex-free
- ✓ Skin-friendly contact layer that meets ISO Standard 10993 for skin irritation and sensitization

Procedure masks meet the following performance criteria:

		Series 1	Series 2
F2219	Sub-micron, particulate filtration efficiency at 0.1 micron %	>95	>98
F1862	Resistance to penetration by synthetic blood, minimum pressure in mm Hg for Pass results	80	120
16 CFR Part 1610	Flame Spread	Class 1	Class 1

Per the FDA's Emergency Use Authorization

Item #	Item Description	Material	Color	Size	Case
88589	Graham Metal Free Procedure Mask	Nonwoven	Blue	OSFM	50/BX; 6BX/CS
86899	Graham Procedure Mask Series 1	Nonwoven	Blue	OSFM	50/BX; 6BX/CS
87529	Graham Procedure Mask Series 2	Nonwoven	Blue	OSFM	50/BX; 6BX/CS

- Not intended to replace the need for FDA-cleared surgical masks or FDA-cleared or authorized respirators.
- Not intended to provide protection against pathogenic biological airborne particulates and are not recommended for use in aerosol generating procedures and any clinical conditions where there is significant risk of infection through inhalation exposure.
- The product has not been FDA cleared or approved.
- The product has been authorized by FDA under an EUA for use in healthcare settings by HCP as PPE to provide a physical barrier to fluids and particulate materials to prevent HCP exposure to respiratory droplets and large particles during surgical mask shortages resulting from the COVID-19 pandemic.
- The product is authorized only for the duration of the declaration that circumstances exist justifying the authorization of the emergency use of medical devices, during the COVID-19 outbreak, under section 564(b)(1) of the Act, 21 U.S.C. § 360bbb-3(b)(1) unless the authorization is terminated or revoked sooner.

Contact Graham Medical for more information on our PPE Product offering!

Disposable, single-use, nonsterile procedure masks are not for operating room use

Graham Medical • 2273 Larsen Road • P.O. Box 19031 • Green Bay, WI 54307-9031
800.558.6765 • info@grahammedical.com • www.grahammedical.com

