

TECHNICAL DATA SHEET

Graham Medical Item #:	54419
Brand:	Graham Medical® PowerFit®
Product Type:	Fitted Sheet
Rating:	NA
Color:	Blue
Material:	Spunbond
Usage:	Disposable
Sterility:	Non-Sterile
Biocompatibility:	NA
Chemotherapy Tested:	NA
Latex-Free:	Yes
Fluid Barrier:	Fluid Resistant
Fastening Type:	Elastic
Size Dimensions*:	33" x 89"
Packaging:	5/BG; 10BG/CS
Country of Origin:	China
Harmonized Code:	6302.32.2040
UNSPSC Code:	42132102



***Note:** To calculate the minimum fitted sheet size needed:

Width = Mattress Width + Twice The Thickness + 2".

Length = Mattress Length + Twice The Thickness + 4".

Graham Medical • 2273 Larsen Road • Green Bay, WI 54303 • 800.558.6765