

## TECHNICAL DATA SHEET

<b>Graham Medical Item #:</b>	54420
<b>Brand:</b>	Graham Medical® PowerFit®
<b>Product Type:</b>	Fitted Sheet
<b>Rating:</b>	NA
<b>Color:</b>	Blue
<b>Material:</b>	Spunbond/Poly
<b>Usage:</b>	Disposable
<b>Sterility:</b>	Non-Sterile
<b>Biocompatibility:</b>	NA
<b>Chemotherapy Tested:</b>	NA
<b>Latex-Free:</b>	Yes
<b>Fluid Barrier:</b>	Fluid Barrier
<b>Fastening Type:</b>	Elastic
<b>Size Dimensions*:</b>	33" x 89"
<b>Packaging:</b>	5/BG; 10BG/CS
<b>Country of Origin:</b>	China
<b>Harmonized Code:</b>	6302.32.2040
<b>UNSPSC Code:</b>	42132102



**\*Note:** To calculate the minimum fitted sheet size needed:

Width = Mattress Width + Twice The Thickness + 2".

Length = Mattress Length + Twice The Thickness + 4".

Graham Medical • 2273 Larsen Road • Green Bay, WI 54303 • 800.558.6765