

TECHNICAL DATA SHEET

Graham Medical Item #:	87475
Brand:	Graham ProDefense™
Product Type:	Anti-Skid Shoe Cover
Rating:	NA
Color:	Blue
Material:	Spunbond
Usage:	Disposable
Sterility:	NonSterile
Chemotherapy Tested:	No
Latex-Free:	Yes
Fluid Barrier:	Fluid-Resistant
Fastening Type:	NA
User:	Adult
Size:	XL
Packaging:	100/BG; 10 BG/CS
Country of Origin:	China
Harmonized Code:	6307.90.9889
UNSPSC Code:	42131609



Graham Medical • 2273 Larsen Road • Green Bay, WI 54303 • 800.558.6765

Revised: 02/08/21