

graham ProDefense™ Protective Gowns

Reduce the risk of cross-contamination and **increase peace of mind** with **disposable** Personal Protective Equipment (PPE).

Poly & Spunbond/Poly



#86792



#86795

- ✓ Long-sleeve, thumb-loop
- ✓ Over-the-head, easy-tear neck
- ✓ Integrated waist tie and open back
- ✓ Poly gown available in two sizes
- ✓ Latex-free

Nonwoven



#86684



#86793

- ✓ Long-sleeve, latex-free elastic wrist
- ✓ Tie neck
- ✓ Waist tie
- ✓ Soft, linen-like material that is fluid-resistant, lightweight and breathable

360° Wrap-Around



#78187



#78627

- ✓ 360° Protection with wrap-around design
- ✓ Latex-free knit cuffs
- ✓ Hook & loop shoulder closure
- ✓ Soft, linen-like material that is fluid-resistant, lightweight and comfortable

Item #	Item Description	Material	Color	Size	Case	Quality	AAMI
86792	Protective Gown - Thumb Loop, OTH Easy-Tear Neck, Waist Tie, Open Back	Poly	Blue	OSFM 42" x 46"	15/BG; 5BG/CS	Premium	Level 2
86900	Protective Gown - Thumb Loop, OTH Easy-Tear Neck, Waist Tie, Open Back	Poly	Blue	XL 46" x 49"	1/BG; 75BG/CS	Premium	Level 2
86795	Protective Gown - Thumb Loop, OTH Easy-Tear Neck, Waist Tie, Open Back	Spunbond/Poly	Blue	OSFM 42" x 46"	50/CS	Premium	Level 2
86793	Isolation Gown - Elastic Wrists, Sewn Neck & Waist Ties	Spunbond	Yellow	M/L 30" x 42"	50/CS	Economy	
86684	Isolation Gown - Elastic Wrists, Sewn Neck & Waist Ties	Nonwoven SMS	Blue	M/L 30" x 42"	50/CS	Standard	
78187	Isolation Gown - 360° Wrap-Around, Knit Cuffs, Hook & Loop Closure	Nonwoven SMS	Yellow	OSFM	25/CS	Standard	
78627	Isolation Gown - 360° Wrap-Around, Knit Cuffs, Hook & Loop Closure	Nonwoven SMS	Dark Green	OSFM	25/CS	Standard	

Contact Graham Medical for more information on our PPE Product offering!

Graham Medical • 2273 Larsen Road • P.O. Box 19031 • Green Bay, WI 54307-9031
800.558.6765 • info@grahammedical.com • www.grahammedical.com

