

TECHNICAL DATA SHEET

Graham Medical Item #:	LabMates® Lab Coats & Jackets
Brand:	Graham Medical
Product Type:	Lab Coats & Jackets
Country of Origin:	China
Harmonized Code:	6210.10.5000
Color:	Blue or White
Material:	Nonwoven SMS
Protection Level:	None - See Note Below *
Chemotherapy Tested:	No
Sterility:	NonSterile
User:	Adult
Usage:	Disposable
Back Style:	Closed
Sleeve Style:	Long Sleeves
Cuff Style:	Knit
Neck Closure:	Knit
Fastening Type:	Snaps
Gender:	Unisex
Size:	Refer to Catalog
Latex-Free:	Yes
UNSPSC Code:	42131607
Packaging:	10/BG; 5BG/CS



* The material used for these products has been tested for water resistance. AATCC42 Impact Penetration: 0.1 gram. AATCC127 Hydrostatic Pressure: 40 cm H₂O.

Graham Medical • 2273 Larsen Road • Green Bay, WI 54303 • 800.558.6765