

graham ProDefense™ MRI Safe Face Mask

Reduce the risk of cross-contamination and **increase peace of mind** with **disposable** Personal Protective Equipment (PPE).

MRI Safe Face Masks Level 1



- ✓ 3-Ply Construction
- ✓ Soft, comfortable ear loops
- ✓ Contains no metal – safe for use in MRI Machines
- ✓ Latex-free
- ✓ Skin-friendly contact layer that meets ISO Standard 10993
- ✓ ASTM Level I Rating



Item Description	Material	Color	Size	Case
Item #88589 Graham ProDefense™ Face Mask MRI Safe Level 1	Nonwoven	Blue	OSFM	50/BX; 6BX/CS



Contact Graham Medical for more information on our PPE Product offering!

Graham Medical • 2273 Larsen Road • P.O. Box 19031 • Green Bay, WI 54307-9031
800.558.6765 • info@grahammedical.com • www.grahammedical.com

